Resin-based Dental Restorative Material

ESTELITE FLOW QUICK

REF PE01

ENGLISH

Carefully read all information, precautions, notes and materials before using ESTELITE FLOW QUICK.

■PRODUCT DESCRIPTION AND GENERAL INFORMATION

- ESTELITE FLOW QUICK is a low viscosity, medium flow, light cured, radiopaque composite resin. Direct placement of ESTELITE FLOW QUICK into a prepared cavity makes handling easy.
- ESTELITE FLOW QUICK adopts Radical-Amplified Photopolymerization initiator technology (RAP technology). Curing time for ESTELITE FLOW QUICK can be shortened by one-third compared with that for our conventional product with a light curing unit having the camphorquinone (CQ) wavelength range (peak: 470nm, spectrum: 400 to 500nm). Please see the table depicting the relationship between curing time and increment depth (■INDICATIONS FOR FILLING AND CURING).
- ESTELITE FLOW QUICK contains 71% by weight (53% by volume) of silica-zirconia filler and silica-titania filler. Every inorganic filler contained in ESTELITE FLOW QUICK is a spherical filler (mean particle size: 300nm, particle size range: 40 to 600nm) that enables excellent gloss retention and wear resistance. The monomer matrix contains Bisphenol A polyethoxy methacrylate (Bis-MPEPP), Triethylene glycol dimethacrylate (TEGDMA), and 1, 6-bis(methacryl-ethyloxycarbonylamino)trimethyl hexane (UDMA).

■SHADE

- A1, A2, A3, A3.5, A4 (Cerv., B1, B2, B3, B4, C1, C2, C3, OA1, OA2, OA3, BW (Bleaching White) and CE (Inc.).
- ESTELITE SIGMA QUICK SHADE GUIDE provided for ESTELITE SIGMA QUICK can also be utilized for ESTELITE FLOW QUICK. Every shade provided for ESTELITE FLOW QUICK is adjusted to that of ESTELITE SIGMA QUICK.
 - Opalescent shades (OA1, OA2, OA3) have adequate opacity designed for blocking out the dark shine that comes through the oral cavity (Class III restoration). The opalescent shades also serve as a dentin layer for the multi-shade layering technique. In addition, the shades OA1, OA2 and BW closely match those of primary teeth. However, they are not designed for masking metal or a darkly stained area such as tetracycline stained teeth.
 - Please check with your local distributor for available shades.

■INDICATIONS

- Direct anterior and posterior restorations (particularly for small/shallow/tunnel shaped cavities)
- As the first layer of direct restorations
- Blocking out cavity undercuts before fabricating indirect restorations
- Repair of porcelain/composite

■CONTRAINDICATIONS

ESTELITE FLOW QUICK contains methacrylic monomers. ESTELITE FLOW QUICK is contraindicated in patients allergic to or hypersensitive to methacrylic and related monomers and should NOT be used in such patients.

■PRECAUTIONS

- DO NOT use ESTELITE FLOW QUICK for any purpose other than those listed in these instructions. Use ESTELITE FLOW QUICK only as directed herein.
- ESTELITE FLOW QUICK is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
- DO NOT use ESTELITE FLOW QUICK if the safety seals are broken or appear to have been tampered with.
- If ESTELITE FLOW QUICK causes an allergic reaction or oversensitivity, discontinue its use immediately.
- The use at all times of examination gloves (plastic, vinyl or latex) is recommended to avoid the possibility of allergic reactions from methacrylic monomers. Certain ingredients may penetrate through examination gloves. If a portion of ESTELITE FLOW QUICK comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- AVOID CONTACT WITH EYES. In case of contact with eyes, immediately flush eyes with water and seek medical attention.
- When using a light-curing unit, protective eye shields, glasses or goggles should be worn.
- In the event that ESTELITE FLOW QUICK comes into contact with the skin, clean the skin immediately with an alcohol soaked cotton swab.
- ESTELITE FLOW QUICK should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- DO NOT leave ESTELITE FLOW QUICK unsupervised within the reach of patients and children to avoid being swallowed by mistake.
- Clean the filling instruments with alcohol after use.

■INDICATIONS FOR FILLING AND CURING

Be sure to light-cure ESTELITE FLOW QUICK extraorally and check the time needed for complete hardening of ESTELITE FLOW QUICK with your light curing unit before performing the bonding procedure. Following table illustrates the relationship between curing time and increment depth.

Relationship between curing time and increment depth:

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■STORAGE


- Store ESTELITE FLOW QUICK at temperatures between 0 - 25°C (32 - 77 °F).
- AVOID direct exposure of ESTELITE FLOW QUICK to light and heat
- DO NOT use ESTELITE FLOW QUICK after the indicated expiration date on the syringe/ package.

■DISPOSAL

Safely dispose of excess ESTELITE FLOW QUICK by extruding unused product from the syringe and light cure before disposal.

■ATTACHING THE SYRINGE TIP (Metal Syringe Tip)

- Hold the nozzle and turn the cap counterclockwise for removal. [1]
- Hold the nozzle and place the Syringe Tip onto it. Turn the Syringe Tip clockwise until it is securely locked. [2]

- Accessories and Syringe Tip are exclusively designed for ESTELITE FLOW QUICK. Avoid using Syringe Tips from other brands with ESTELITE FLOW QUICK.
- Syringe Tips are also sold separately.
- Use Gloves when attaching the Syringe Tip.

■CLINICAL PROCEDURE

- Cleaning**
Teeth should be cleaned with pumice and water to remove surface stains and plaque.

- Shade Selection**
Select the appropriate shade of ESTELITE FLOW QUICK.

- Isolation**
A rubber dam is the preferred method of isolation.

- Cavity Preparation**
Prepare the cavity. In case where no cavity preparation has been made (caries-free cervical defects), clean the tooth surface with a rubber cup and a fluoride-free cleaning paste. Rinse thoroughly with water.

- Pulp Protection**
Glass ionomer lining or calcium hydroxide should be applied if the cavity is in close proximity to the pulp. DO NOT USE EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of bonding agent and/or ESTELITE FLOW QUICK.

- Bonding system**
Apply light-cured bonding system according to its manufacturer's instructions.
 - Do not use self-cured bonding systems. We advise that light-cured bonding systems should be used with ESTELITE FLOW QUICK since other types of bonding systems such as self-cured or dual-cured bonding systems may not provide effective adhesion when used with light-cured composite resins. If you want to use self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.
 - In case of porcelain/composite repairs be sure to pretreat the surface of porcelain/composite with a silane coupling reagent BEFORE performing the bonding procedure.

- Filling and contouring**
Cavities can be directly filled from the syringe, or indirectly filled with an instrument after extruding the paste on the mixing pad. Increments should not exceed the indicated depth at a time (please refer to the aforementioned table).
 - DO NOT apply unnecessary force to the syringes, immediately after removal from the refrigerator.
 - DO NOT mix ESTELITE FLOW QUICK with other brands of composite resin to avoid incomplete cure or air bubble contamination. In order to avoid air bubble contamination, DO NOT mix ESTELITE FLOW QUICK with another shade of the paste.
 - After extruding the paste, remove and discard the Syringe Tip, wipe the nozzle with gauze and attach the cap immediately.

- Curing**
Light-cure each increment for at least the indicated time (please see aforementioned table).
 - If other brands of composite resins are layered over the cured composite, follow their instructions.

- Finishing**
Shape and polish the restoration. For finishing, use fine finishing diamond points. To erase the demarcation between the composite and the natural tooth rotate a 12 fluted carbide bur slowly along the margins without water. Use finishing metal strips or vinyl strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tool, such as PoGoTM (DENTSPLY/Caulk), Sof-LexTM (3M ESPE), Identoflex® HiLuster Dia Polishers (Kerr), or D-FINETM Hybrid Diamond (COSMEDENT, Inc.).
 - *PoGo, Sof-Lex, Identoflex and D-FINE are not registered trademarks of Tokuyama Dental Corporation.

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the dental professional to ensure the product is suitable for application before use. Specifications are subjected to change without notification.

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