



# VITA MFT

Multi Functional Teeth

Mould chart acrylic teeth




Date of issue 12.08




**VITA**

Exclusively available in the  
VITA SYSTEM 3D-MASTER shades



Ovoid



Triangular

**040**



**T41\***



**044**




**T43**



**049\***



**T46**



**T49**



\* Available 06/2009




**VITA MFT – The new dimension in aesthetics**

The natural-looking anatomy of VITA MFT anteriors, with approximal contacts in the area of the incisal third being slightly placed towards the palate, enables a trouble-free and appearingly organic alignment of the teeth.

The anatomically correct labial curve leads to a physiologically perfect set-up of the teeth where discretely broadened tooth necks support an age-appropriate papilla design.

## Lower anteriors



Rectangular

**R42**



**R45**



**R47\***




Square

**S47**



**S50\***



**L33**



**L34**



**L37**



**L37LN\***



**L40\***



The mould **L37LN** (Long Neck) was specifically developed for simulating particular periodontal characteristics which, for example, occur in partial prosthetics.

\* Available 06/2009

**PU29 – PL29**



**PU31 – PL31**




**PU33 – PL33**



**easy-centric® – the new, effortless way for setting up dentures**

The occlusal surface design of VITA MFT is based on the innovative principle called easy-centric which supports a professional and economical set-up procedure. Moreover, it leads to increased productivity in a quality well adapted to the patient's needs.

No.		No.	mm	
	mm	mm	No.	mm
<b>O40</b>	40.0	<b>L33</b> 32.6	<b>PU29</b> 28.8	<b>PL29</b> 29.6
<b>O44</b>	43.9	<b>L34</b> 34.0	<b>PU29</b> 28.8	<b>PL29</b> 29.6
<b>O49</b>	49.2	<b>L40</b> 39.5	<b>PU33</b> 33.0	<b>PL33</b> 33.6
<b>T41</b>	40.5	<b>L33</b> 32.6	<b>PU29</b> 28.8	<b>PL29</b> 29.6
<b>T43</b>	43.3	<b>L37</b> 36.8	<b>PU31</b> 30.6	<b>PL31</b> 31.1
<b>T46</b>	46.1	<b>L37</b> 36.8	<b>PU31</b> 30.6	<b>PL31</b> 31.1
<b>T49</b>	48.8	<b>L40</b> 39.5	<b>PU33</b> 33.0	<b>PL33</b> 33.6
<b>R42</b>	42.1	<b>L34</b> 34.0	<b>PU29</b> 28.8	<b>PL29</b> 29.6
<b>R45</b>	44.5	<b>L37</b> 36.8	<b>PU31</b> 30.6	<b>PL31</b> 31.1
<b>R47</b>	46.6	<b>L37</b> 36.8	<b>PU33</b> 33.0	<b>PL33</b> 33.6
<b>S47</b>	46.8	<b>L40</b> 39.5	<b>PU33</b> 33.0	<b>PL33</b> 33.6
<b>S50</b>	49.7	<b>L40</b> 39.5	<b>PU33</b> 33.0	<b>PL33</b> 33.6











#### VITA MFT – Anteriors and posteriors

The whole range of indications in the pleasantly small VITA MFT product line. A reduced and harmonic range of shades/moulds makes for increased economic efficiency.

The multifunctional VITA MFT tooth is suited for:


- Full dentures
- Implant prosthetics
- Telescope and cone constructions
- Attachment Dentures
- Model cast prostheses
- Interim solutions

The shade range of VITA MFT covers the 8 most common VITA SYSTEM 3D-MASTER shades:

VITA SYSTEM 3D-MASTER							
							
1M1	2L1.5	2M1	2M2	3L2.5	3M2	3R2.5	4M2
(~ B1)*	(~ B2)*	(~ A1)*	(~ A2)*	(~ B3)*	(~ A3)*	(~ A3.5)*	(~ A4)*

\* Best Match VITAPAN classical

With the unique VITA SYSTEM 3D-MASTER, all natural tooth shades can be systematically determined and perfectly reproduced.



Please note: Our products should be used according to the working instructions. We cannot be held liable for damages resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with porcelains and equipment from other manufacturers which are not compatible or not authorized for use with our product. Furthermore, our liability for the correctness of this information is independent of the legal ground and, in as far as legally permissible, is limited to the invoiced value of the goods supplied excluding turnover tax. In particular, as far as legally permissible, we do not assume any liability for profit loss, for indirect damages, for consequential damages or for claims of third parties against the purchaser. Claims for damages based on fault liability (culpa in contrahendo, breach of contract, unlawful acts, etc.) can only be made in the case of intent or gross negligence. Date of issue of these directions for use: 09.08

VITA Zahnfabrik is certified according to the Medical Device Directive; the following products bear the mark  0124: VITA MFT

US 5498157 A · AU 659964 B2 · EP 0591958 B1

# VITA

VITA Zahnfabrik H. Rauter GmbH & Co.KG  
Postfach 1338 · D-79704 Bad Säckingen · Germany  
Tel. +49(0)7761/562-0 · Fax +49(0)7761/562-299  
Hotline: Tel. +49(0)7761/562-222 · Fax +49(0)7761/562-446  
[www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com) · [info@vita-zahnfabrik.com](mailto:info@vita-zahnfabrik.com)