



Elastomerics



Global Experts in Orthodontics

03

GC Orthodontics
GLOBAL EXPERTS IN ORTHODONTICS

Elastomerics Summary

Elastomeric color chains	143
Orchestra Intra and extra-oral Elastics	144
Japanese elastomeric chain "GC Ortho Chain"	145
GC Ties elastomeric ligatures	145
Separators	146
Hollow or solid elastomeric wires	146
Archwire protection sleeves	147



Elastomeric color chains

Available in a wide range of colors that remain stable in time, our color chains long or short – are sold in 4.5 m spools.



	Color	Short	Long
■	Silver	88-1001-0000	88-3001-0000
■	Pearl	88-1002-0000	
□	Clear	88-1003-0000	88-3003-0000
■	Blue	88-1004-0000	88-3004-0000
■	Red	88-1005-0000	
■	Green	88-1006-0000	88-3006-0000
■	Black	88-1007-0000	88-3007-0000
□	White	88-1008-0000	
■	Pink	88-1009-0000	
■	Metallic Purple	88-1010-0000	
■	Metallic Blue	88-1011-0000	
■	Metallic Red	88-1012-0000	
■	Lavender	88-1013-0000	88-3013-0000
■	Fuchsia	88-1014-0000	
■	Yellow	88-1015-0000	
■	Orange	88-1016-0000	88-3016-0000
■	Navy Blue	88-1017-0000	
■	Light Blue	88-1018-0000	88-3018-0000
■	Maroon	88-1019-0000	
■	Grey	88-1020-0000	88-3020-0000
■	Glow Blue	88-1021-0000	
■	Metallic Gold	88-1022-0000	88-3022-0000



Orchestra Intra and extra-oral Elastics

Our range of intra and extra oral elastics is made from the highest quality medical-grade latex.

The stated force corresponds to the force the elastic exerts when stretched up to 3 times its original diameter. Bags of 100 intra-oral elastics, 25 bags per box. Bags of 50 extra-oral elastics, 25 bags per box. P/N in white in the table are latex-free elastics.

				EXTRA-ORAL
Size \ Force	2.5 oz (70,88 g)	4.0 oz (113,40 g)	6.0 oz (170,10 g)	8,5 oz (240,98g)
3 mm (1/8")	85-0321-0000 Harp	85-0331-0000 Xylophone	85-0341-0000 Tuba	–
4 mm (3/16")	85-0421-0000 Tambourine	85-0431-0000 Piano	85-0441-0000 Trumpet	85-0461-0000 Trombone
4 mm (3/16") (latex free)	–	85-0430-0000 Jazz (latex free)	85-0440-0000 Rock-n-Roll (latex free)	–
6 mm (1/4")	85-0621-0000 Oboe	85-0631-0000 Flute	85-0641-0000 Baritone	85-0661-0000 French Horn
6 mm (1/4") (latex free)	–	85-0630-0000 Hip-Hop (latex free)	85-0640-0000 Classical (latex free)	–
8 mm (5/16")	85-0821-0000 Cello	85-0831-0000 Clarinet	85-0841-0000 Saxophone	85-0861-0000 Guitar
10 mm (3/8")	85-1021-0000 Cymbals	85-1031-0000 Tympani	–	85-1061-0000 Violin



Japanese elastomeric chain "GC Ortho Chain"

Made from high-density and high-purity materials, this chain is precision-cut and intervals between links are strictly identical. All along its action, it exerts light and continuous forces, essential requirements for optimal tooth movements. Available in grey or clear. In 2m spools.

Description	P/N
GC Ortho Chain Clear Short	86-1010-0000
GC Ortho Chain Clear Medium	86-1020-0000
GC Ortho Chain Clear Long	86-1030-0000
GC Ortho Chain Grey Short	86-2010-0000
GC Ortho Chain Grey Medium	86-2020-0000
GC Ortho Chain Grey Long	86-2030-0000




GC Ties elastomeric ligatures

Latex-free, GC Ties ligatures will allow an upper and lower archwires change with a single stick.

The quality of the polyurethane enables these ligatures to recover 90 % of their initial shape after being stretched.

Available in bags of 50 sticks, 20 ligatures per stick.

Outer diameter: .120"

 Silver	87-1201-0000	 Lavender	87-1213-0000
 Pearl	87-1202-0000	 Fuchsia	87-1214-0000
 Clear	87-1203-0000	 Yellow	87-1215-0000
 Blue	87-1204-0000	 Orange	87-1216-0000
 Red	87-1205-0000	 Navy Blue	87-1217-0000
 Green	87-1206-0000	 Light Blue	87-1218-0000
 Black	87-1207-0000	 Maroon	87-1219-0000
 White	87-1208-0000	 Grey	87-1220-0000
 Pink	87-1209-0000	 Glow Blue	87-1221-0000
 Metallic Purple	87-1210-0000	 Golden	87-1222-0000
 Metallic Blue	87-1211-0000	 Glow Clear	87-1223-0000
 Metallic Red	87-1212-0000	 Tooth	87-1224-0000



Separators

Latex-free. Very strong, with a memory effect, these separators are ideal for creating spaces between molars. Sold by 1000.

Description	P/N
radiopaque separators	89-1602-0000

Hollow or solid elastomeric wires



Our hollow or solid elastomeric wires are available in two sections: .025" and .030". 15m spools.

Hollow wires:	P/N
.025" Grey	89-4325-0000
.025" Clear	89-3325-0000
.030" Clear	89-3335-0000

Solid wires:	P/N
.025" Grey	89-4125-0000
.025" Clear	89-3125-0000
.030" Grey	89-4135-0000
.030" Clear	89-3135-0000

Solid fiber with treated surface	P/N
.025" Clear	89-3222-0000
.025" Grey	89-4222-0000
.030" Grey	89-4232-0000



Archwire protection sleeves

3m spools

For archwires up to .018"	P/N
Grey	89-4311-0000
Clear	89-3311-0000

For archwires up to .021" x .027"	P/N
Grey	89-4321-0000
Clear	89-3321-0000

For .045" archwires	P/N
Grey	89-4341-0000

